

# Trilogy100

## Portable homecare ventilator

Trilogy100 is a life-support ventilator designed to provide invasive and noninvasive home ventilation for both adults and children weighing at least 11 pounds.

### Easy to use

Trilogy100 provides easy-to-read, easy-to-navigate screens with clear, concise directions.

- The device can switch between detailed displays for clinicians and simplified screens for patients and caregivers.
- The intuitive design allows for quick access to device settings and patient information.

### Portable

Trilogy100's portable design provides patients with increased mobility and makes it easy for clinicians and caregivers to move the device.

- Size: 6.6" x 11.2" x 9.3"
- Includes a 3-hour internal battery and a 3-hour detachable battery that can be quickly exchanged with spare batteries.
- Is lightweight (11 pounds) and has fewer connections than a standard circuit.

### Versatile

Trilogy100 is both a pressure and volume control ventilator.

- Features Respironics' proven BiPAP technology with digital Auto-Trak for true leak compensation.
- Allows clinicians to set dual prescriptions—one for daytime and a second for when the patient is asleep.
- Interchangeable ports allow clinicians to choose best available circuit and patient interface.
- Can accommodate a mask, mouthpiece, or tracheostomy.

### DirectView software

To assist clinicians, DirectView patient data management software is available for use with Trilogy100. DirectView streamlines reporting to help clinicians identify adjustments necessary to maintain effective treatment and improve patient care. DirectView enables quick access to full patient therapy information including waveforms, trends, usage patterns, and summary statistics for broader views of patient therapy.

